



**SURGICAL
INNOVATIONS**

HAND & FOOT

Innovative Simplicity



Z-Medical

We offer | Innovative Simplicity



Z-Staples

International patented **EU/US/CHN**



Z-SnapOff Screws

(Compression Screws flat head)

International patented **EU/US/CHN**



Z-Screws cannulated / headless

(Double Compression / Round head compression / Neutralisation)



Z-Screwdriver Technology

International patented **EU/US/CHN**



- › Various USPs/product provide advantages in sales & marketing
- › “Ready to use” sterile implants with 7 years shelf-life
- › Full traceability by UDI
- › Reusable instruments available



**SURGICAL
INNOVATIONS**



Z-SnapOff Screws


Innovative Simplicity



Innovative Simplicity | Z-SnapOff Screws

Features

Snap Off Screws
Flat head
Compression



∅ 2mm

Indication

- › Internal bone fixation in case of bone fractures, bone reconstruction, fusions or osteotomies in ankles, feet and hands.
- › Especially fixation of small bone fragments, weil osteotomy, mono-cortical fixation, osteotomies and fractures fixation in the foot and hand.


Z-Medical


Innovative Simplicity | Z-SnapOff Screws



SURGICAL
INNOVATIONS

Z-SnapOff Screws with Navigator or Z-Speed Drive system





- › Trocar tip enables multidirectional application
- › Self drilling
- › Self tapping
- › Thread profile for optimized compression

- › Innovative and patented SnapOff principle to pre-fixation
- › Increase side-stability secures angled insertion
- › Flash-free breakage of the application shaft

- › TORX drive to endposition

Innovative Simplicity | Z-SnapOff Screws (Flat Head)



C01 015 Ø2 11mm



C01 017 Ø2 13mm



C01 016 Ø2 12mm



C01 018 Ø2 14mm

- › Implants sterile packed
- › Titanium alloy



STERILE

+



Reusable instruments for all types of screws available!

=

READY TO USE

Z-Medical

Innovative Simplicity | Z-SnapOff Screws



SURGICAL
INNOVATIONS

Reusable instrumentation available!




Z-Screw insertion “SnapOff”
(power unit)



Z-Screw Navigator
(screw positioning and length measurement)



Z-Screw Driver
(direct screw insertion)



Screw Driver

Screw Navigator

Z-Speed Drive System „Manual”

- › Screw positioning
- › Screw length determination
- › Screw insertion



Screw insertion
with power unit

Z-Speed Drive System „power unit”

- › Screw insertion “Snap Off”



Features

- › Patented Z-SnapOff technology
(directional stability against early snap-off)
- › Speed Drive System for screw insertion by power unit
- › Screw Navigator for screw positioning available
- › Sterile packed implants available
- › Reusable instruments available



**SURGICAL
INNOVATIONS**

Z-Medical® GmbH + Co. KG
Gänsäcker 38, 78532 Tuttlingen
Germany

fon: +49 (0) 7462 94 55 40

fax: +49 (0) 7462 94 55 49

email: info@z-medical.de

web: www.z-medical.de

Follow us:



Z-Medical® GmbH + Co. KG